



THE NEW GENERATION FOR LIVE CELL SHIPPERS

CELLBOX
Ground CD



The Cellbox Ground CD makes it possible to perform very long distance transport by road, rail and sea. It was designed to use disposable CO₂ cartridges, allowing for quick and easy refilling on-the-go. The user-programmable CO₂ and temperature values allow you to optimize the conditions in the incubation chamber for your precious biological material.

Benefits of live cell transport



- ➔ No pH-shift or use of metabolically active cryoprotectants
- ➔ No altered protein expression; no mechanic or osmotic stress
- ➔ Intelligently maintains internal environment for up to 24 hours
- ➔ High cell viability and accurate continuous data logging
- ➔ No freezing, no thawing - Cells arrive ready-to-use



CO₂ SUPPLY
Provides regulated CO₂ environment



TEMPERATURE
Maintains temperature between 28-38°C



APP
Data logging and export via Cellbox App



RECHARGEABLE
Li-ION battery with external 100 - 230V power supply

CONTACT

+49 (0) 40 226 316 4-10

enquiries@cellbox-solutions.com

www.cellbox-solutions.com

Englische Planke 8 20459 Hamburg | An der Bahn 5 22844 Norderstedt Germany